

# Urine Control With Biochemical Diagnostics

## Detectabuse™ AU/NZ MIX For Hand Held Devices



### Top Reasons to use Detectabuse™ Controls:

- Human urine-based matrix simulates patient samples
- Recommended by the majority of hand held device manufacturers
- 12 drug constituents within one control vial
- Stable for up to 3 years frozen, or refrigerated (unopened)
- Ready to use in 5 mL vial sizes
- Use Detectabuse™ controls to train staff in the accurate use of hand held devices
- Use Detectabuse™ controls to evaluate the performance of various devices on the market
- Each Detectabuse™ control contains: THC, Benzoylcgonine (COC), PCP, Morphine (OPI), Methamphetamine, Amphetamine, Secobarbital (BAR), Oxazepam (BZO), Methadone, Methaqualone, Propoxyphene, Nortriptyline (TCA).

### The Detectabuse Advantage:

Detectabuse AU/NZ MIX controls are formulated to work with all known hand held single use drug-screening devices. We also offer a complete line of controls for immunoassay and GC/MS

### NEW! AUSTRALIAN-NEW ZEALAND CUTOFF LEVELS

TARGET VALUES	-50%	Cutoff*	+50%
Delta-9-THC-COOH	25	50	75
Benzoylcgonine	150	300	450
Phencyclidine (PCP)	12.5	25	37.5
Morphine	150	300	450
d-Methamphetamine	150	300	450
d-Amphetamine	150	300	450
Secobarbital	150	300	450
Oxazepam	100	200	300
Methadone	150	300	450
Methaqualone	150	300	450
Propoxyphene	150	300	450
Nortriptyline	500	1000	1500

jump to Pkg.  
Insert

### Detectabuse AU/NZ MIX Multi-Constituent Urine Control

19700000	Negative	5 mL	\$8.40 USD
19480010	Cutoff -50%	5 mL	\$8.40 USD
19480021	Cutoff +50%	5 mL	\$8.40 USD



Contact: Colleen Gang, Product Manager, Liquid Controls • 180 Heartland Blvd, Edgewood, NY 11717  
 • E-mail: cgang@biochemicaldiagnostics.com • Phone: 1-800-223-4835 • Fax: 1-631-595-9204