



Biochemical Diagnostics, Inc.

180 Heartland Blvd., Edgewood, NY 11717 • Phone (631) 595-9200 Fax (631) 595-9204

MATERIAL SAFETY DATA SHEET

GENERAL INFORMATION

TRADE NAME: Sodium Chloride

Catalog Number: 3507100-8, 3507105-9,
3506213-1, 3506214-9, 3507000-1,
Part of Kits: 3506303-0, 3507020-6,
3506112-6, 3506261-1, 3506094-4

Chemical Name: Sodium Chloride

CAS #: N/A

Formula: NaCl

INGREDIENTS

Components: Sodium Chloride

Threshold Limit Values: N/A

PHYSICAL DATA

Boiling Point: 1413 C/2575 F

Vapor Pressure: <1 mmHg

Vapor Density: N/A

Specific Gravity: 2.160

Percent Volatile: 0%

Evaporation Rate: N/A

Solubility H₂O: Appreciable >10%

pH: N/A

Appearance: White granules or solid

Odor: Odorless

Melting Point: 801 C/1473 F

FIRE AND EXPLOSION HAZARD DATA

Flammable Limits in air % by volume:

Lel: N/A

Uel: N/A

Flash Point: N/A

Auto Ignition Temp: N/A

Extinguishing Media: Use extinguishing media appropriate for surrounding fire

Special Fire Fighting Procedures: Wear proper protective equipment w/SCBA

Unusual Fire and Explosion Hazards: None identified

REACTIVITY DATA

Stability: Stable

Conditions to Avoid: None documented

Incompatibility: Strong oxidizing agent, lithium, interhalogens (ex: bromine, trifluoride), iron, cement.

Hazardous Polymerization: Will not occur

Hazardous Decomposition: Hydrogen chloride

ENVIRONMENTAL PROTECTION PROCEDURES

Spill Response: Wear suitable clothing. Absorb with vermiculite.

Waste Disposal Method: Dispose in accordance with state, local and federal regulations.

SPECIAL PROTECTION INFORMATION

Eye Protection: Safety goggles

Skin Protection: Gloves, Lab Coat

Respiratory Protection: None Required

Ventilation: Use adequate general or local exhaust.

SPECIAL PRECAUTIONS

Insert Precautions:

Keep container tightly closed.

Contaminated Equipment Precautions: Wash thoroughly with soap and water.

Hygienic Handling and Storage: Wash thoroughly with soap and water.

EMERGENCY AND FIRST AID PROCEDURES:

Ingestion: May be harmful if swallowed in large quantity.

Inhalation: Remove to fresh air.

Skin Contact: Wash immediately with soap and water.

Eye Contact: Rinse with copious amounts of water.

Health:	1	Slight
Flammability:	0	None
Reactivity:	1	Slight
Contact:	1	Slight

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, we make no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independent judgment in determining its appropriateness for a particular purpose. Accordingly, we will not be responsible for damages of any kind resulting from the use of or reliance upon such information. No representations, or warranties, either express or implied of merchantability fitness for a particular purpose or of any other nature made hereunder with respect to the information set forth herein or to the product to which the information refers.