

# Biochemical Diagnostics Detectabuse Liquid Urine Controls for 2008 SAMHSA CUTOFFS

- Human urine-based matrix simulates patient samples
- Up to 12 drug constituents within one control
- Stable for up to 3 years frozen, or refrigerated (unopened)
- Series X-H in ready to use 20 mL vial size.
- Confirm 2 in ready to use 50 mL vial size.
- Each Series X-H control contains: THC, Benzoylcgonine, PCP, Morphine, 6-MAM, d-Methamphetamine, MDMA, Secobarbital, Oxazepam, Methadone, Methaqualone, Propoxyphene
- Each Confirm 2 control contains: THC, Benzoylcgonine (COC), PCP, Codeine, Morphine (OPI), 6-MAM, d-Amphetamine d-Methamphetamine, MDMA, MDA, MDEA



[SERIES-X-H  
Pkg. Insert](#)



[CONFIRM-2  
Pkg. Insert](#)

The DETECTABUSE™ Advantage:

Detectabuse™ Series X-H & Confirm 2 controls are designed to monitor the performance of qualitative or quantitative procedures for the detection of drugs in urine. We also offer a complete line of controls for other immunoassay, hand held single use devices, GC/MS and custom formulations.

NEW! DETECTABUSE™ LIQUID CONTROL URINE IMMUNOASSAY SERIES X-H AND GC/MS CONFIRM 2			
Catalog #	Series X-H (High Opiate, 2000 ng/mL Cutoff)	Vial Size	List Price
19601052	Series X-H Cutoff -25%	20 mL	\$29.80
19601056	Series X-H Cutoff +25%	20 mL	\$29.80
19601908	Series X-H 3X Cutoff (special Order Only)	20 mL	\$50.40
GC/MS Confirm 2 (High Opiate, 2000 ng/mL Cutoff)			
19582001	Confirm 2 Cutoff -60%	50 mL	\$74.50
19582501	Confirm 2 Cutoff +25%	50 mL	\$74.50

[Go to pricelist](#)



Contact: Colleen Gang, Product Manager, Liquid Controls • 180 Heartland Blvd, Edgewood, NY 11717  
 • E-mail: [cgang@biochemicaldiagnostics.com](mailto:cgang@biochemicaldiagnostics.com) • Phone: 1-800-223-4835 • Fax: 1-631-595-9204  
 BCD West – Contact: Maria Hunt, Manager • Phone: 1-800-223-8803 • E-mail: [mhunt@biochemicaldiagnostics.com](mailto:mhunt@biochemicaldiagnostics.com)